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Centrifuges with nitrogen purging

Application and instrument configuration

Nitrogen flows through the centrifuge chamber and reduces its oxygen concentration. Suitability for the sample medium, solvents, vessels and operating conditions must be assessed for the specific instrument configuration. Nitrogen purging alone does not authorise centrifugation of flammable or explosive substances and does not constitute explosion protection certification of the centrifuge.

Purging and start release

Before purging, the motor and the relevant electrical consumers are switched off. Nitrogen enters from below; purge gas is discharged at the top. Overpressure monitoring and an oxygen probe monitor the chamber atmosphere. If the wrong gas supply is connected, such as compressed air instead of nitrogen, monitoring does not release the centrifuge for starting.

Monitoring and shutdown

The overpressure monitoring system and oxygen probe have separate power supplies. If either system reports a fault, a safety relay de-energises the relevant current-carrying components. The rotor then coasts unbraked to a standstill. This shutdown function does not replace assessment of the specific application.

Power failure

Bypass purging can maintain nitrogen flow during a mains power failure. This requires the corresponding instrument configuration and an available nitrogen supply.

Components and configuration scope

The underlying technical description from November 2024 identifies the overpressure system, oxygen probe and safety relay as SIL 2 components. This statement concerns the described components; it does not establish overall SIL 2 or explosion protection certification of the complete centrifuge. Components, monitoring limits, gas connections and protective functions must be confirmed for the specific configuration.

Information for selection and quotation

Please specify the centrifuge model, sample medium and solvents, vessels and closures, required rotor, speed or RCF, temperature and available gas supply. Scope of supply and interfaces are defined in the quotation; function illustrations do not establish the scope of supply.

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