


Accessories and Spare Parts

UV Protection



Safety Advice:

Always wear face/eye protection
when working with UV light!

Also protect your skin from UV
light.



Cat. No.	Description
29 59 298	UV Protection Shield with hinges for Economy Transilluminators UVT-14 and UVT-22
29 70 098	UV Protection Shield with hinges for Economy Transilluminators UVT-20 and UVT-28
29 93 100	UV Protection Glasses
29 93 000	UV Full Face Protection

Accessories Blue Light

Cat. No.	Description
29 93 150	Orange Glasses



All products built by Herolab carry the CE sign. Herolab is certified according to DIN EN ISO 9001:2008 and DIN EN ISO 13485.
All prices without VAT, ex works. – Herolab reserves the right to change technical specifications and prices without prior notice.



Replacement Tubes



for Transilluminators

For Model		Energy	Wave Length	No. of Tubes per Unit	Cat. No.
UVT-14 UVT-22	S	6 Watt	254 nm	4	29 84 800
	M	6 Watt	312 (302) nm	4	29 84 900
	L	6 Watt	365 nm	4	29 85 000
UVT-20	S	8 Watt	254 nm	6	29 84 300
	M	8 Watt	312 (302) nm	6	29 84 400
	L	8 Watt	365 nm	6	29 84 500
UVT-20	S	8 Watt	254 nm	6	29 84 300
	M	8 Watt	312 (302) nm	6	29 84 400
	L	8 Watt	365 nm	6	29 84 500
Combi units	SML	6 Watt	White Light	je 3	29 84 600
UVT-28	S	8 Watt	254 nm	8	29 84 300
	M	8 Watt	312 (302) nm	8	29 84 400
	L	8 Watt	365 nm	8	29 84 500
UVT-40	S	15 Watt	254 nm	6	29 84 000
	M	15 Watt	312 (302) nm	6	29 84 100
	L	15 Watt	365 nm	6	29 84 200
White light units	W	15 Watt	White Light	2	29 84 700
UVT-14 BE		18 Watt	Blue Light	1	29 84 650
UVT-20 BE		18 Watt	Blue Light	2	29 84 650
UVT-28 BE		18 Watt	Blue Light	3	29 84 650

All products built by Herolab carry the CE sign. Herolab is certified according to DIN EN ISO 9001:2008 and DIN EN ISO 13485.
All prices without VAT, ex works. – Herolab reserves the right to change technical specifications and prices without prior notice.



Replacement Tubes for UV Handlamps

For Model		Power	Wavelength	No. of Tubes per Unit	Cat. No.
UV-4	S	4 Watt	254 nm	1	29 83 300
	L	4 Watt	365 nm	1	29 83 500
	BLB	4 Watt	365 nm	1	29 83 508
UV-6 UV-12	S	6 Watt	254 nm	1 resp. 2	29 84 800
	M	6 Watt	312 (302) nm	1 resp. 2	29 84 900
	L	6 Watt	365 nm	1 resp. 2	29 85 000
	BLB	6 Watt	365 nm	1 resp. 2	29 85 008
UV-8 UV-16	S	8 Watt	254 nm	1 resp. 2	29 84 300
	M	8 Watt	312 (302) nm	1 resp. 2	29 84 400
	L	8 Watt	365 nm	1 resp. 2	29 84 500
	BLB	8 Watt	365 nm	1 resp. 2	29 84 508
UV-15 UV-30	S	15 Watt	254 nm	1 resp. 2	29 84 000
	M	15 Watt	312 (302) nm	1 resp. 2	29 84 100
	L	15 Watt	365 nm	1 resp. 2	29 84 200
	BLB	15 Watt	365 nm	1 resp. 2	29 84 208

Other accessories and spare parts on request.

All products built by Herolab carry the CE sign. Herolab is certified according to DIN EN ISO 9001:2008 and DIN EN ISO 13485.
All prices without VAT, ex works. – Herolab reserves the right to change technical specifications and prices without prior notice.

