

➤ **New:** Centrifuge Bottles with wide mouth - 70 mm

Material: PPCO

PPCO-Bottles 70 mm neck	Volume	Dimensions Tested up to	Cat. No. Bottle with Cover	Cat. No. Bottle w./o. Cover	Cat. No. Cover
 blue cap	750 ml	Ø 100 x 143 mm For Fixed angle rotor 8,000 rpm / 12,042 x g	25 37 30	23 37 30	25 46 80
	950 ml	Ø 100 x 167 mm For Fixed angle rotor 8,000 rpm / 12,042 x g	25 37 40	23 37 40	25 46 80
	1,000 ml	Ø 100 x 181 mm For Fixed angle rotor 8,000 rpm / 12,042 x g	25 37 50	23 37 50	25 46 80
 white cap	1,000 ml	Ø 98 x 190 mm For Fixed angle rotor 8,000 rpm / 12,042 x g	25 37 80	23 37 80	25 46 73
 black cap	750 ml	Ø 100 x 130 mm For Swing-out rotor 4,100 rpm / 5,394 x g	25 37 35	23 37 30	25 46 81
	950 ml	Ø 100 x 154 mm For Swing-out rotor 4,100 rpm / 5,394 x g	25 37 45	23 37 40	25 46 81
	1,000 ml	Ø 100 x 168 mm For Swing-out rotor 4,100 rpm / 5,394 x g	25 37 55	23 37 50	25 46 81

All products built by Herolab carry the CE sign. Herolab is certified after DIN EN ISO 9001:2008 and DIN EN ISO 13485.
All prices without VAT, ex works. – Herolab reserves the right to change technical specifications and prices without prior notice.