



Centrifuge HiCen with pressurized enclosure according DIN EN 60079-2



- 1) Centrifuge Chamber
- 2) Rotor
- 3) Rotor-Lid
- 4) Centrifuge Tube/Bottle
- 5) Tube Closure
- 6) Pressure valve at the outlet
- 7) Pressure measurement
- 8) Oxygen measurement
- 9) Over-temperature safety switch (Bimetal)
- 10) Temperatur sensor
- 11) ID-Sensor
(Rotor-in, Rotor-Type, Overspeed monitoring)
- 12) Inlet proportional valve
- 13) Safety valve
(Throttle check valve, normally open)
- 14) Outlet
- 15) Sealing chamber lid
- 16) Gasket to motor

E. Knorr 1.3.2018